



Band and Loop Space Maintainer

The following is about Band and Loop Space Maintainer related, I hope to help you better understand Band and Loop Space Maintainer.

China Band and Loop Space Maintainer outsourcing dental lab ,manufacturer and Dental implants dental center Called wmdentallab



We have now a skilled, performance group to offer excellent support for our consumer. We usually follow the tenet of customer-oriented, details-focused for factory Outlets for Shenzhen Wanmei Dental Laboratory Co.,LTD.

Band and Loop Space Maintainer, Dental Unit, Denture, Orthodontic Products, Dental Products, From China Dental Lab, Any requires from you will be paid with our best attention!

factory Outlets for China Dental, Dental Product, top quality, timely delivery and your satisfaction are guaranteed. We welcome all inquiries and comments. If you are interested in any of our items or have an OEM order to fulfill, remember to feel free to contact us now. Working with us will save you money and time.

The Band and Loop appliance is a unilateral space maintainer that is used to prevent the drifting of adjacent teeth in edentulous areas after exfoliation or extraction. Often used in mixed dentition cases to maintain space to allow for eruption of a permanent tooth. May be designed for unilateral or bilateral treatment.

When a patient loses a primary molar thereby leaving a space in the dentition, the first molar tends to drift forward into that space causing crowding and blocking the proper eruption of the bicuspid(s).

A Space Maintainer (also know as a Band & Loop) uses a wire loop soldered to a molar band and extending anteriorly where it contacts the distal edge of the tooth on the mesial side of the space to prevent this drift. The appliance is left in place until the patient's bicuspid(s) erupt and fill the space.